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Water based Surface Disinfection

FIELD OF APPLICATION

VERIFORTE® is a rapid and effective water-based disinfectant for use on surfaces. VERIFORTE® is pH neutral, alcohol-free, aldehyde-free, fragrance-free and non-flammable.

INSTRUCTIONS FOR USE

VERIFORTE® is a ready to use, fluid, clear solution. Surfaces and objects are simultaneously and thoroughly moistened with 50ml/m² using where appropriate a disposable and leaving to dry. Rinsing with water is not necessary. The product evaporates leaving no residue.

HIGHLY EFFECTIVE

VERIFORTE® is tested according to European Standards EN 14885, EN 13727, EN 13624, EN 13697, EN 13704, EN 1276 and EN 14476 and is listed by DGHM/VAH. VERIFORTE® can be used both as a wipe and spray disinfectant. The stated exposure times below, demonstrate the full effectiveness of Veriforte against bacteria, viruses, spores and yeast, even when, mechanical help is not employed.

EXPOSURE TIMES undiluted without mechanical help:

Bacterial Effect	surface precleaned 5min*	
	MRSA.....5 min	E.coli.....1 min
	E.hirae.....1 min	P.aeruginosa.....1 min
Yeasticidal Effect	precleaned surface	
	C. albicans.....5 min	
	*shortest possible exposure time according to VAH – Applied Hygiene Network	
Sporocidal Effect	precleaned surface	
	Bac. Subtilis.....5 min	
Virucidal Effect	Murine Norovirus	low burden 1 min / high burden 1 min
	Adenovirus	low burden 1 min / high burden 1 min
	Poliovirus	low burden 1 min / high burden 30 min
Virucidal against	enveloped viruses.....5 min according to the DVV and RKI guidelines	
(incl. HIV, HBV und HCV)		

SPECIFICATIONS Color/Form: transparent/liquid pH value: 6.5---7.5

COMPOSITION >99.5% Water; <0.4% NaCl (Sodium Chloride); <0.05% resp. 0.5g/L NaOCl (Sodium Hypochlorite) / HOCl (Hypochlorous acid)

REGISTRATION
 N – 54547 for Product type 2 (disinfectants)
 N – 54587 for Product type 4 (Food and feed sector)
 N – 54588 for Product type 5 (drinking water)

HAZARDS IDENTIFICATION

The product is non-hazardous. Normal precautions should however be observed when handling chemicals: Note the lot number and expiration date printed on the packaging.

SAFETY

In normal use, VERIFORTE® is not dangerous for people. Avoid eye contact. Irritation may occur. After eye contact, rinse immediately with clear water for several minutes (remove contact lenses). If you experience lengthy complaints related to the product, please consult a medical specialist.





STORAGE

Store VERIFORTE® in a light-protected and cool place.
 Keep container closed tightly. Note the lot number and expiration date printed on the packaging.

DISPOSAL

With a sufficient dilution Veriforte can be disposed via a standard wastewater sewage system. Please observe local and country-specific waste and wastewater regulations. After emptying and thorough rinsing with water, the plastic container can be recycled or can be disposed of in an appropriate way.

SYMBOLS

Use biocides safely. Always read the label and product information before use.

Manufacturer: P.G.F. Industry Solutions GmbH | Katzmoosstraße 26a | 5161 Elixhausen | Austria | office@veriforte.com | +43(662) 846540

The information in this product data sheet corresponds to our current level of knowledge. They are the result of our previous experience and scientific reviews. However, we cannot individually examine and influence different circumstances in practical application. The quality of the work therefore depends on the professional processing by the user. The user is obligated to always check the compatibility of the material to be treated with the product in an inconspicuous place or to obtain technical advice. Liability for incomplete or incorrect information in our data sheet exists only in case of gross negligence or intent. Any claims arising from the Product Liability Act remain unaffected. With this product data sheet, all previous product information will lose their validity. Document Name: DI-7-205-3, 26.05.2020